Meta-analysis shows patient wearable sensor reduces incidence of hospital acquired pressure injuries in critically ill patients

Leo Nherera, Annemari Cooley, Gary Delhougne

Background and Aims

- Each year, more than 2.5 million people in the United States develop pressure ulcers , which are the second most common but deadliest Hospital-Acquired Condition (HAC), yet the majority of them are preventable.
- Several best practices have been shown to be effective in reducing the occurrence of pressure ulcers, most notably regular and frequent repositioning of immobile patients around the clock.
- However, current methodologies such as paper clocks or musical chimes to prompt on-time repositioning have not been demonstrated to produce long-lasting, sustainable compliance to turn protocols and therefore sustained outcomes.
- A wearable patient sensor system (patient sensor) has been shown to provide an objective recording of patient turning practices.
- The purpose of this study was to assess the clinical effectiveness of patient sensor in the prevention of hospital acquired pressure injuries (HAPIs) in acutely ill patients.

Materials and methods

Literature search

- A targeted non-systematic literature review was performed to identify published studies and conference abstracts
- We searched conference websites and used contacts from the marketing department to identify studies

Inclusion Criteria

- Comparative studies- comparing the patient sensor versus standard of care
- Full publication or conference abstract
- Report on incidence of HAPIs

Meta-analysis

• Meta-analyses were performed; a fixed effect model was used when heterogeneity was not significant ($I^2 < 50\%$) and a random effects model when it was significant ($I^2 \ge 50\%$) to estimate the overall effect of the patient sensor in reducing the HAPIs

Results

- Literature search identified 8 published studies
- 1 randomised controlled trial (RCT) and 7 conference abstracts
- A total of 34,711 patients were included in the studies
- 19,136 patients used the patient sensor.
- The majority of the patients were over 75 years of age.
- The patient sensor reduced incidence of HAPIs by 70% Risk ratio (RR) 0.30 (95% Confidence interval 0.21- 0.44) p<0.00001.

Discussion & Conclusions

- Over 19,000 patients provided data on the effectiveness of the patient sensor in reducing HAPIs.
- The meta-analysis evidence suggests that the sensor reduces HAPIs by 70%.
- Evidence from one RCT reported the reduction at 73%.
- More controlled studies are needed to validate these results.

A meta-analysis of 19,136 patients demonstrated that the patient sensor* reduced the incidence of hospital acquired pressure injuries by 70% (p < 0.00001)

Meta-analysis results, Incidence of HAPI, patient sensor versus standard care

	Patient sensor		Standard of care			Risk Ratio
Study or Subgroup	Event	Total	Event	Total	Weight	M-H, Random, 95% Cl
Rosini 2020	4	875	24	809	8.4%	0.15 (0.05, 0.44)
Freshcorn 2018	18	2172	52	2172	17%	0.35 (0.20, 0.59)
Hendricken 2019	3	597	99	597	6.2%	0.33 (0.09, 1.23)
Larson 2018	21	3475	14	555	14.1%	0.24 (0.12, 0.47)
McManus 2017	27	9823	38	9258	17.9%	0.67 (0.41, 1.10)
Ohnstat 2018	6	1270	19	1270	10.1%	0.32 (0.13, 0.79)
Pickham 2018	5	671	15	555	8.9%	0.28 (0.10, 0.75)
Rogers 2020	16	253	85	273	17.6%	0.20 (0.12, 0.34)
Total (95% CI)		19136		15489	100%	0.30 (0.21, 0.44)

100 256 Total events Heterogeneity: Tau^z=0.14; Chi^z= 14.50, df= 7 (P=0.04); I^z=52% Test for overall effect: Z= 6.30 (P<0.00001)

> References 1. Agency for Healthcare Research and Quality: Preventing Pressure Ulcers in Hospitals. https://www.ahrq.gov/ patient-safety/settings/hospital/resource/pressureulcer/tool/pu1.html retrieved on September 30, 2019

> 2. Agency for Healthcare Research and Quality: AHRQ National Scorecard on Hospital-Acquired Conditions Updated Baseline Rates and Preliminary Results 2014–2016. June 2018. Retrieved from https://www.ahrq.gov/ sites/default/files/wysiwyg/professionals/quality-patient-safety/pfp/natlhacratereport-rebaselining 2014-2016_0.pdf on September 30, 2019

* Leaf Patient Monitoring system (Leaf Healthcare, Inc. Pleasanton, CA)

SmithAephew



©2020 Smith + Nephew, Inc. ♦Trademark of Smith + Nephew. LFTE1-22453-0120